

Jazz Ultra

With you every step of the way



Keeler



Small in size, big on features

Go anywhere

Weighing in at just 89g - less than half the weight of an apple - and measuring just 130mm, the Jazz Ultra is the perfect fit for your pocket or handbag.



Lighter

Whether you are walking the hospital wards or the corridors of your surgery, the 'go anywhere' Jazz Ultra, with its sleek ergonomic design, will be with you every step of the way.

Illuminating

With the introduction of LED illumination and a colour temperature of natural daylight, the Jazz Ultra lights the way in ophthalmoscopy and otoscopy. The newly designed optical chassis provides greater optical quality and illumination accuracy.



Brighter



Supplied with 3 interchangeable colour rings



Green for life

One AA battery

Only requiring one AA battery, the Jazz Ultra's LEDs have been designed to be extremely energy efficient.

This means you can have the freedom to use the Jazz Ultra all day without worrying about battery life.



Greener

Long life LEDs

Designed to last more than 20,000 hours, the LED in the Jazz Ultra will last more than 15 years*. Making the Jazz Ultra a once in a lifetime purchase.



Long lasting

* Based on 4 hours use per day.



Ophthalmoscope

- Brilliant LED illumination and trusted Keeler optics provide the clearest, crispest views of the retina
- Lens powers +/- 1, 2, 3, 4, 6, 8, 10, 15 and 20 Dioptres
- Smooth aperture wheel operation for precise selection
- Protects your spectacles

Apertures and Filters



- 1 On/Off switch
- 2 Dioptre display
- 3 Dioptre adjustment wheel
- 4 Aperture wheel
- 5 Interchangeable colour rings
- 6 Spectacle protection
- 7 Handle/battery compartment





*Instruments shown
at actual size*



*Only one AA
battery is required
to power Jazz Ultra
instruments*

Otoscope

- Pure white uniform illumination shows colours as they really are for more efficient diagnosis plus a longer life battery and no more replacement of bulbs
- Sealing precision optics with a triple magnification for optimum views every time
- Can be used for performance of pneumatic tests (with connector and insufflator)
- Specula available in 5 sizes: 2mm, 2.5mm, 3mm, 4mm and 5mm

- ① On/Off switch
- ② 3x magnifying pivoting lens
- ③ Pneumatic test connection
- ④ Fibre optics
- ⑤ Interchangeable colour rings
- ⑥ Ear speculum
- ⑦ Handle/battery compartment

Part Numbers and Accessories



**Jazz Ultra Pocket
Ophthalmoscope
with Pouch**

Part Number **1102-P-1113**



**Jazz Ultra Pocket
Diagnostic Set
in Double Pouch**

Part Number **1702-P-1109**



**Jazz Ultra Pocket
Otoscope with Pouch**

Part Number **1514-P-1187**



Ophthalmoscope LED

Part Number **1011-P-7122**



Otoscope LED

Part Number **1015-P-7074**



Spare Coloured Rings

Part Numbers

EP39-37051 Black Ring

EP39-37342 Pink Ring

EP39-37350 Blue/Grey Ring



Disposable Specula

Part Numbers

1514-P-7056 2mm Pack of 10

1514-P-7044 2.5mm Pack of 10

1514-P-7052 3mm Pack of 10

1514-P-7060 4mm Pack of 10

1514-P-7079 5mm Pack of 10

1514-P-7087 2mm Pack of 100

1514-P-7095 2.5mm Pack of 100

1514-P-7108 3mm Pack of 100

1514-P-7116 4mm Pack of 100

1514-P-7124 5mm Pack of 100



Go anywhere



Distributor